

### Cum sa utilizezi corect dispozitivul Turbuhaler



1. Desfaceti si inlaturati sigiliul de protectie



2. Tineti dispozitivul in pozitie verticala. Rotiti mineral cat este posibil, intr-o directie si apoi in sens inver



3. Expirati. Nu expirati prin piesa bucala.



4. Introduceti piesa bucala intre dinti, strangeti buzele si inspirati adanc si cu putere prin piesa bucala



5. Scoateti dispozitivul din gura inainte de a expira



6. Puneti la loc capacul prin insurubare stransa